Anti-M monoclonal, mouse

Deutsch

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Customer Service:
service@bag-diagnostics.com
Tel.: +49 (0) 6404/925 - 125

Anti-M monoclonal, mouse

For Tube Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-M reagent is produced from cell culture supernatants of hybridoma-cell lines. The cells are secreting antibodies of IgG-type which reacts specific with the corresponding blood group antigen. The reagent is used for In-Vitro-Diagnostic, to determine qualitative, whether human red blood cells possess or lack the corresponding blood group antigen M.

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the appropriate antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed reagent contains antibodies of the following clone:

Anti-M monoclonal, mouse clone: M-11H2

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagent is prepared from supernatants of cell cultures. Regardless, as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagent to indicated expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. A weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required.
Take and use reagent directly from the vials.

PROCEDURE

Not provided material, additionally needed

- Tube Centrifugation Method:
1. tubes (10 x 75 mm or 12 x 75 mm)
 2. microliter pipette
 3. timer
 4. centrifuge
 5. isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure**Tube Centrifugation Method**

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the reagent in a marked tube, subsequently add 100 µL of appropriate cell suspension.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL reagent.
3. Mix Erythrocytes-/Reagent-mixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifuge tube for 1 min at 1.000 rpm (approximately 180-270 x g).
6. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method:

- Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.
Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. Due to variability of antigen expression, reactivity of the reagent against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

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| | | | |
|---|------------------------------|---|-----------------------------|
|  REF | Product Code |  LOT | Lot |
|  | Store from - to |  | Expiration Date |
|  IVD | In-Vitro Diagnostic |  | EU CE symbol |
|  | Manufacturer as to 98/79/EG |  | Observe instruction for use |
|  UDI | Unique Device Identification |  | Distributor |

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BAG Diagnostics GmbH