AHG serum blend

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AHG serum blend

Coombs serum, polyspecific, polyclonal, rabbit Anti-IgG / monoclonal, mouse Anti-C3 clone: Bric-8

For direct and indirect Coombs-Test
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Polyspecific Coombs serum for blood group serology contains antibodies against Immunoglobulins and complement factors. The present AHG serum blend contains polyclonal reagent from the plasma of immunized rabbits (Anti-IgG) and monoclonal reagent from cell culture supernatants from mouse hybridoma cell line Bric-8 (anti-C3). The form of preparation takes care, that no reaction with uncoated human erythrocytes takes place.

For secure addition of reagent green colour is added.
The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method use with this reagent is based on the principle of agglutination. Antigens present on human erythrocytes are directly or indirectly detected and loaded with the corresponding antibody of the IgG-Type. The antibodies of the AHG serum blend (secondary antibody), reacts with the IgG immunoglobulins and / or C3 complement factors (C3d and C3b) that are bound to the surface of the red blood cells.

REAGENTS

The listed polyspecific Coombs serum blend is available in following formulation:

AHG serum blend Coombs serum, polyspecific, polyclonal, rabbit Anti-IgG / monoclonal, mouse Anti-C3 clone: Bric-8

The reagent contains <0.1% (w/v) sodium azide as preservative.
Beside the parts active antibody and animal serum, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagent is prepared from animal plasma and supernatants of cell cultures.
As biological products they should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.
On disposal, flush with large quantities of water.
For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C.
May be at room temperature while in use.
In principle, store and use the reagent to indicated expiry date only.

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagent.
- Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
- Strength of positive reactions also depends on age of used blood
- Centrifugation outside the specified speed range may lead to false results.
- The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of this reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

SAMPLE PREPARATION

- Blood sample should be collected by approved medical procedure.
- Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required.
Take and use reagent directly from the vials.

PROCEDURE

Not provided material, additionally needed

- Tube Centrifugation Method:
- tubes (10 x 75 mm or 12 x 75 mm)
 - microliter pipette
 - timer
 - incubator
 - centrifuge
 - isotonic saline (0,85 - 0,9% sodium chloride)
 - Anti-Human-Globulin-Reagent (Coombs-Serum / AHG- Serum)

Test procedure

Direct Anti-Globulin-Test (Coombs-Test)

Tube Centrifugation Method

- Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
- At first put 100 µL of the AHG serum blend in a marked tube, subsequently add 100 µL of appropriate cell suspension.
Alternative one drop = approximately 50 µL cell suspension can be added to two drops = approximately 100 µL AHG serum blend.
- Mix Erythrocytes- / Reagent-mixture well by slightly shaking.
- Incubate tube at room temperature for 5-15 min.
- Centrifuge tube for 1 min at 1.000 rpm (approximately 180-270 x g).
- Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
- Document the result.

Indirect Anti-Globulin-Test (Coombs-Test)

Tube Centrifugation Method

- Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
- At first put 100 µL of reagent in a marked tube, subsequently add 100 µL of appropriate cell suspension.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL reagent.
- Mix Erythrocytes- / Reagent-mixture well by slightly shaking.
- Incubate tube in an incubator at +37 °C for 30 min.
- Wash red cells 3 times with (cold) isotonic saline.
- Subsequently add 100 µL AHG serum blend (Coombs-Serum) to the tube, release the cells from the bottom of the tube by slightly shaking and mix with the Coombs-Serum.
- Centrifuge tube for 1 min at 1.000 rpm (approximately 180-270 x g).
- Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
- Document the result.

INTERPRETATION OF RESULTS

"Slightly shake up" of the erythrocyte sediment at Tube Centrifugation Method.

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and / or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted samples should not be used.
- Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may give weaker reactivity compared to control cells.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.²
- Insufficient washing of erythrocytes may cause false results.
- Contamination of reagent with human protein may inactivate the antibodies and cause false negative results
- As described in the literature, samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵

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REF	Product Code	LOT	Lot
	Store from - to		Expiration Date
IVD	In-Vitro Diagnostic	CE	EU CE symbol
	Manufacturer as to 98/79/EG		Observe instruction for use
UDI	Unique Device Identification		Distributor

REF

- 69255** AHG serum blend Coombs serum, polyspecific, polyclonal, rabbit Anti-IgG / monoclonal, mouse Anti-C3 clone: Bric-8 10 ml
- 69256** AHG serum blend Coombs serum, polyspecific, polyclonal, rabbit Anti-IgG / monoclonal, mouse Anti-C3 clone: Bric-8 10 x 10 ml

730-13-1709 Version 009 / 20. January 2022

BAG
DIAGNOSTICS

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering:
T. +49 (0) 6404 / 925 - 450
F. +49 (0) 6404 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com