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Anti-C, Anti-c, Anti-E, Anti-e monoclonal, human IgM

English

For Slide-, Tube-, Card- and Microplate-Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-C, Anti-c, Anti-E and Anti-e reagents are produced from cell culture supernatants of heterohybridoma - cell lines. The cells are secreting an antibody of IgM-type that reacts specific with the respective corresponding blood group antigen. The antibody in each case is human protein. The reagents are used for In-Vitro-Diagnostic, to determine qualitative whether human red blood cells possess or lack the corresponding blood group antigens C, c, E, and e.

The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with the reagents mentioned above are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed monoclonal, human IgM reagents are available in two different variations. The differences depends only in the antibody-producing clones:

Anti-C clones: MS-24, P3x25513G8	Anti-C clone: MS-273
Anti-c clone: MS-33	Anti-c clone: MS-35
Anti-E clones: MS-258, 906	Anti-E clones: MS-12, MS-260
Anti-e clones: MS-16, MS-21, MS-63	Anti-e clones: MS-62, MS-69

These reagents contain 0.1% (w/v) sodium azide as preservative. Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents are prepared from supernatants of cell cultures. As biological products they should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagents to indicated expiry date only.

REMARKS

- With each testing positive and negative controls should be performed
- Inappropriate storage impairs efficacy of the reagents.
- Weak turbidity of the reagents mentioned above does not affected its reactivity. Bacteria and chemical contamination of the products should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
- Strength of positive reactions also depends on age of used blood.
- Centrifugation outside the specified speed range may lead to false results.
With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of the reagents mentioned above all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“ 1 in its current version.
- The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

- Blood sample should be collected by approved medical procedure.
- Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the samples.
If a delay in testing occurs, samples should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required.
Take and use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed

Slide Method:	glass slide; Pasteur pipette; mixing stick; timer;
Tube Centrifugation Method:	test tubes (10 x 75 mm or 12 x 75 mm); microliter pipette; centrifuge; timer; isotonic saline (0,85 - 0,9% sodium chloride)
Card Method:	Cards: - BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“ - Grifols „DG Gel Neutral“; - Ortho BioVue® System Reverse; test tubes; microliter pipette; timer; centrifuge; corresponding card centrifuge; card specific diluent; isotonic saline (0,85 - 0,9% sodium chloride)
Microplate Method:	microplate with 96 U-wells; microliter pipette; centrifuge; timer; microplate shaker; microplate centrifuge; isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

- Use erythrocyte sediment only.
- Place one drop (approximately 50 µL) of appropriate reagent on a marked glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
- Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
- Document the result.

Tube Centrifugation Method

- Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
- At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
- Mix Erythrocytes- / Reagent mixture well by slightly shaking.
- Incubate tube at room temperature for 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 - 1.000 x g).
- Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
- Document the result.

Card Method (manual method / valid for the cards:

- BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“
 - Grifols „DG Gel Neutral“
- Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
 - Add 50 µL of appropriate cell suspension to a marked micro tube.
 - Add 25 µL of appropriate reagent to the micro tube.
 - Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force after latest 30 minutes.
 - Check macroscopically micro tube for agglutination within 30 minutes.
 - Document the result.

Ortho BioVue® System

- Prepare 3-5 % suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
- Add 40 µL of appropriate reagent to a marked reaction chamber.
- Add 10 µL of appropriate cell suspension to the reaction chamber.
- Centrifuge the cassette after latest 15 minutes in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
- The test results should be read macroscopically for agglutination directly after the end of centrifugation.
- Document the result.

Microplate Method

- Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
- Add 50 µL of appropriate reagent to a marked well.
- Add 50 µL of appropriate cell suspension to the well.
- To mix both shake microplate for 30 seconds on microplate-shaker with medium speed.
- Centrifuge microplate in appropriate microplate-centrifuge for 30 seconds at 400 x g.
- Shake microplate for 30 seconds on a microplate shaker with medium speed.
- Evaluated macroscopically for agglutination directly after shaking.
- Document the result.
- Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
- Repeat steps 5 to 8 after incubation at room temperature.

INTERPRETATION OF RESULTS

"Slightly rotating / shaking"

at Slide Method / at Tube Centrifugation Method / Microplate Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card according to the instruction for use of the corresponding card.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted samples should not be used.
- With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
- Due to variability of antigen expression, reactivity of the reagents mentioned above against certain phenotypes may give weaker reactivity compared to control cells.
- Some red blood cells may have developed Rh variants that may not be recognized with monoclonal antibodies of a clone. Therefore the tests should / have to be carried out with at least two different clones / cloning mixtures as required in some countries. Deviating test results require further testing.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens².
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in card method. These cells react positive without test serum too.
- Pay attention to all statements to limitations in the instruction for use of the cards.
- Negative test results with anti-e (clones: MS-62, MS-69) in the microplate method has to be confirmed with another technique, as false negative results may occur.
- Do not use blood of newborn on microplates.
- With blood of newborn anti-c (clone: MS-33) may possess anti-I/i cold agglutinin activity and therefore show false positive results.

BAG
DIAGNOSTICS

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering:
T. +49 (0) 6404 / 925 - 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products

Product	Positive	Negative	False Positive	False Negative
Anti-C clones: MS-24, P3x25513G8	649	305	0	0
Anti-C clone: MS-273	649	305	0	0
Anti-c clone: MS-33	777	177	0	0
Anti-c clone: MS-35	777	177	0	0
Anti-E clones: MS-258, 906	249	705	0	0
Anti-E clones: MS-12, MS-260	249	705	0	0
Anti-e clones: MS-16, MS-21, MS-63	938	16	0	0
Anti-e clones: MS-62, MS-69	1138	63	0	0*

The calculated values found for:

Product	Sensitivity	Specificity
Anti-C clones: MS-24, P3x25513G8	100%	100%
Anti-C clone: MS-273	100%	100%
Anti-c clone: MS-33	100%	100%
Anti-c clone: MS-35	100%	100%
Anti-E clones: MS-258, 906	100%	100%
Anti-E clones: MS-12, MS-260	100%	100%
Anti-e clones: MS-16, MS-21, MS-63	100%	100%
Anti-e clones: MS-62, MS-69	100%*	100%

* for Anti-e (clones MS-62, MS-69) in microplate method 42 samples were tested false negative, so in this method sensitivity was 97,12% only

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REF Product Code	LOT Batch
 Store from - to	 Expiration Date
IVD In-Vitro Diagnostic	CE EU CE symbol
 Manufacturer according to 98/79/EG	 Consult instructions for use
UDI Unique Device Identification	 Distributor

730-13-0520 Version 020 / 15. January 2022

CE 0483

REF

69214	Anti-C monoclonal, human IgM clones: MS-24, P3X25513G8	10 ml
69215	Anti-C monoclonal, human IgM clones: MS-24, P3X25513G8	10 x 10 ml
69216	Anti-C monoclonal, human IgM clone: MS-273	10 ml
69217	Anti-C monoclonal, human IgM clone: MS-273	10 x 10 ml
69220	Anti-c monoclonal, human IgM clone: MS-33	10 ml
69221	Anti-c monoclonal, human IgM clones: MS-33	10 x 10 ml
69222	Anti-E monoclonal, human IgM clone: MS-35	10 ml
69223	Anti-c monoclonal, human IgM clone: MS-35	10 x 10 ml
69224	Anti-E monoclonal, human IgM clones: MS-258, 906	10 ml
69225	Anti-E monoclonal, human IgM clones: MS-258, 906	10 x 10 ml
69226	Anti-E monoclonal, human IgM clones: MS-12, MS-260	10 ml
69227	Anti-E monoclonal, human IgM clones: MS-12, MS-260	10 x 10 ml
69228	Anti-e monoclonal, human IgM clones: MS-16, MS-21, MS-63	10 ml
69229	Anti-e monoclonal, human IgM clones: MS-16, MS-21, MS-63	10 x 10 ml
69230	Anti-e monoclonal, human IgM clones: MS-62, MS-69	10 ml
69231	Anti-e monoclonal, human IgM clones: MS-62, MS-69	10 x 10 ml