Anti-Jk^a (Ant-JK1) (monoklonal, human IgM)

Anti-Jk^b (Ant-JK2) (monoklonal, human IgM)

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Anti-Jk^a (Ant-JK1) (monoklonal, human IgM)

Anti-Jk^b (Ant-JK2) (monoklonal, human IgM)

For Tube-Test
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-Jk^a and -Jk^b reagents are produced from cell culture supernatants of heterohybridoma-cell lines. The cells are secreting antibodies of IgM-type that reacts specific with the corresponding antigen. The antibody is a human protein. The reagents are used to determine whether red blood cells possess or lack the corresponding blood group antigens Jk^a and Jk^b. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The procedure used with these reagents is based on the principle of agglutination. Normal human erythrocytes, possessing one of these antigens, will agglutinate in the presence of the specific antibody directed toward the antigen.

PRODUCT NAME

The listed reagents contain antibodies of the following cell clones:

Anti-Jk^a (monoclonal, human IgM, clone MS-15)

Anti-Jk^b (monoclonal, human IgM, clone MS-8)

Both reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody the reagents contain sodium chloride, macromolecules and bovine albumin.

WARNING

These reagents were prepared from supernatants of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. Store and use the reagents only to declared expiry date in principle.

REMARKS

- Strength of positive reactions also depends on age of used blood.
- With each testing positive and negative controls should be executed.
- Inappropriate storage impairs efficacy of the reagents.
- Centrifugation different from appointed relative centrifugal force will lead to faulty results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

PROCEDURE

Not provided material additionally needed

at Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver approximately 100 µL, 50 µL
- Centrifuge
- Isotonic saline; 0,85 - 0,9% sodium chloride

Test procedure

Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate reagent to each tube.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube.
- Mix well by slightly shaking.
- Incubate tube at room temperature (15 to 30 °C) for 15 min.
- Centrifugation of tube for 1 min at 2.000 rpm (approximately 800-1.000 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

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LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may produce incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecific.
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

PERFORMANCE CHARACTERISTICS

A performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	Negative	False Positive	False Negative
Anti- Jk ^a	193	170	0	0
Anti- Jk ^b	209	170	0	0

The calculated values found for:

Product	Sensitivity	Specificity
Anti- Jk ^a	100%	100%
Anti- Jk ^b	100%	100%

CE 0483

730-13-1108 Version 008 / 20. Januar 2022

BAG
DIAGNOSTICS

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering:
T. +49 (0) 6404 / 925 - 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com

ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany