# Anti-P<sub>1</sub> monoclonal, mouse

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**Anti-P<sub>1</sub> monoclonal, mouse**

For Slide- and Tube Method  
FOR IN VITRO DIAGNOSTIC USE ONLY

**INTENDED USE**

Monoclonal agglutinating Anti-P<sub>1</sub> reagent is produced from cell culture supernatants of a hybridoma-cell line. The cells are secreting antibodies of IgM type, which reacts specific with the corresponding blood group antigen. The reagent is used for In-Vitro-Diagnostic, to determine qualitative, whether human red blood cells possess or lack the corresponding blood group antigen P<sub>1</sub>.

The reagent is intended to be used by qualified and technical personnel only.

**PRINCIPLE OF PROCEDURE**

The test method used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the appropriate antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

**REAGENTS**

The listed reagent contains antibodies of the following clone:

Anti-P<sub>1</sub> monoclonal, mouse clone: 650

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

**WARNING**

The reagent is prepared from supernatants of cell cultures. Nevertheless, as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

**STORAGE REQUIREMENT**

Store opened and unopened products at +2 to +8°C.

May be at room temperature while in use.

In principle, store and use the reagent to indicated expiry date only.

**REMARKS**

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. A weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagent all effective national laws, directives and guidelines in its current valid version have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup>.

**SAMPLE PREPARATION**

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.  
If a delay in testing occurs, sample should be stored at +2 to +8°C.  
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.  
Blood bag / Donor Blood can be tested until the expiry date.

**REAGENT PREPARATION**

There is no preparation of the reagent required.

Take and use reagent directly from the vials.

**PROCEDURE**

Not provided material, additionally needed

Slide Method:

1. glass slide
2. pasteur pipette
3. mixing stick
4. timer

Tube Centrifugation Method:

1. tubes (10 x 75 mm or 12 x 75 mm)
2. microliter pipette
3. timer
4. centrifuge
5. isotonic saline (0,85 - 0,9% sodium chloride)

**Test procedure****Slide Method**

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of the reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

**Tube Centrifugation Method**

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the reagent in a marked tube, subsequently add 100 µL of appropriate cell suspension.  
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL reagent.
3. Mix Erythrocytes- / Reagentmixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifuge tube for 1 min at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

**INTERPRETATION OF RESULTS**

"slightly rotating" at Slide Method  
"slightly shaking" at Tube Centrifugation Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

**LIMITATIONS OF THE PROCEDURE**

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
6. Due to variability of antigen expression, reactivity of the reagent against certain phenotypes may give weaker reactivity compared to control cells.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.<sup>2</sup>
8. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

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 REF	Product Code	 LOT	Lot
	Store from - to		Expiration Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer as to 98/79/EG		Observe instruction for use
 UDI	Unique Device Identification		Distributor

**REF**

69238 Anti-P<sub>1</sub> monoclonal, mouse clone: 650 5 ml

730-13-2707 Version 007 / 20. January 2022

**BAG**  
DIAGNOSTICS

**BAG Diagnostics GmbH**  
Amtsgerichtsstr. 1-5  
35423 Lich/Germany  
**A BAG Group company**

**T.** +49 (0) 6404/925-100  
**F.** +49 (0) 6404/925-460  
**M.** info@bag-diagnostics.com  
**W.** www.bag-diagnostics.com

**Ordering:**  
**T.** +49 (0) 6404 / 925 - 450  
**F.** +49 (0) 6404 / 925 - 460  
**M.** order@bag-diagnostics.com

**Customer Service:**  
**T.** +49 (0) 6404 / 925 - 125  
**F.** +49 (0) 6404 / 925 - 421  
**M.** service@bag-diagnostics.com