## Anti-Fy<sup>a</sup> Coombs-reaktiv, polyklonal, human

## Anti-Fy<sup>b</sup> Coombs-reaktiv, polyklonal, human

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## Anti-Fy<sup>a</sup> Coombs-reaktive, polyklonal, human

## Anti-Fy<sup>b</sup> Coombs-reaktive, polyklonal, human

For indirect Coombs-test  
FOR IN VITRO DIAGNOSTIC USE ONLY

### INTRODUCTION

Polyclonal Coombs-reaktive Anti-Fy<sup>a</sup> und -Fy<sup>b</sup> - reagents are produced from human plasma that contains a specific antibody of IgG-type, which reacts exclusively with the corresponding antigen. The reagents are used to determine whether red blood cells possess or lack the corresponding blood group antigens Fy<sup>a</sup> and Fy<sup>b</sup>. The reagents are intended to be used by qualified and technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedure used with these reagents is based on the principle of agglutination. Normal human erythrocytes, possessing one of these antigens, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

### PRODUCT NAME

The listed reagents are available in a formulation, with differences in specificity for the antibody only:

Anti-Fy<sup>a</sup> Coombs-reaktive, polyklonal, human

Anti-Fy<sup>b</sup> Coombs-reaktive, polyklonal, human

Both reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody and human serum the reagents contain sodium chloride, macromolecules and bovine albumin.

### WARNING

These reagents were prepared from human plasma. The raw materials for these products are tested for HBsAg, HIV 1/2- and HCV antibodies and found to be negative. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above reagents should be handled with proper care.

### STORAGE REQUIREMENT

Store unopened and opened products tightly sealed at +2 to +8°C. May be at room temperature (+15 to +30°C) while in use. Store and use the reagents only to declared expiry date in principle.

### REMARKS

- Strength of positive reactions also depends on age of used blood.
- With each testing positive and negative controls should be executed.
- Weak turbidity of the reagent does not affect its reactivity.
- Inappropriate storage impairs efficacy of the reagents.
- Centrifugation highly different from appointed relative centrifugal force will lead to faulty results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of these products on such automated systems.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the current „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

### REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

### PROCEDURE

Not provided material additionally needed

#### at Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver approximately 100 µL, 50 µL
- Centrifuge
- Incubator
- Isotonic saline; 0,85 - 0,9% sodium chloride
- Anti-human-globulin-serum (Coombs-serum)

#### Test procedure

##### Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate reagent to each tube.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube.
- Mix well by slightly shaking.
- Incubate tube at +37 °C for 30 min.
- Wash red cells 3 times with (cold) isotonic saline.
- Add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to each tube.
- Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

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## LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may produce incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions with these reagents.
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
6. The described procedures for testing are for usage with Coombs-serum produced by Antitoxin GmbH. On principle Coombs-sera of other manufacturers can be used, but laboratories have to follow the procedures of this manufacturer. Also laboratories must follow approved validation procedures to demonstrate compatibility of these products.

## PERFORMANCE CHARACTERISTICS

A performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	Negative	False Positive	False Negative
Anti- Fy <sup>a</sup>	179	224	0	0
Anti- Fy <sup>b</sup>	201	247	0	0

The calculated values found for:

Product	Sensitivity	Specificity
Anti- Fy <sup>a</sup>	100%	100%
Anti- Fy <sup>b</sup>	100%	100%

## LITERATUR

Peter D. Issitt, David J. Anstee Applied Blood Group Serology, fourth edition, Montgomery Scientific Publications 1998

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BAG Diagnostics GmbH