NUR ZUR IN-VITRO-DIAGNOSTIK

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REF

- 02.184-02** Anti-Jk^a for card method Coombs-reactive, polyclonal, human 2 ml
- 02.184-03** Anti-Jk^a for card method Coombs-reactive, polyclonal, human 3 ml
- 02.184-05** Anti-Jk^a for card method Coombs-reactive, polyclonal, human 5 ml
- 02.185-02** Anti-Jk^b for card method Coombs-reactive, polyclonal, human 2 ml
- 02.185-03** Anti-Jk^b for card method Coombs-reactive, polyclonal, human 3 ml
- 02.185-05** Anti-Jk^b for card method Coombs-reactive, polyclonal, human 5 ml

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FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Polyclonal Coombs-reactive Anti-Jk^a and Anti-Jk^b reagents are produced from human plasma that contains antibody of IgG-type, which reacts specific with the corresponding blood group antigen. The reagents are used to in-vitro determine qualitative whether human red blood cells possess or lack the corresponding blood group antigens Jk^a and Jk^b. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with these reagents are based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

REAGENTS

The listed reagents are available in following formulation:

Anti-Jk^a for card method Coombs-reactive, polyclonal, human

Anti-Jk^b for card method Coombs-reactive, polyclonal, human

The reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody and human serum, the reagents contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagents are prepared from human plasma. The raw material for these reagents are tested for HBsAg, HIV 1/2- and HCV-antibodies and found to be negative. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagents to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagents should no longer be used this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the appropriate card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test methods identified below are for manual testing only and for the devices listed in this insert. When using other automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on such automated systems.
7. For usage of the reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“ in its actual form.
8. The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood samples should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagents.
If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required.
Take and use the reagents directly from the vials.

PROCEDURE

Not included but necessary materials:

Card Method manuell:

1. Cards:
 - Diagnostic Grifols „DG Gel Coombs“
 - Ortho BioVue® System AHG polyspecific cassettes or Ortho BioVue® System Anti-Humanglobulin Anti-IgG cassettes (BIO-RAD (DiaMed) cassette: only applicable with the test serum Anti-Jk^a)
 - BIO-RAD (DiaMed) ID-Card "Liss/Coombs"
2. Microliter pipette
3. Tubes
4. Isotonic saline (0,85 - 0,9% sodium chloride)
5. Centrifuge
6. Card Centrifuge
7. Card incubator
8. Card specific diluent

Card Method automatic:

Vision® Analyzer

Test procedure

Card Method (manual method / valid for the cards:

- Diagnostic Grifols „DG Gel Coombs“)
- (BIO-RAD (DiaMed) cassette: only applicable with the test serum Anti-Jk^a)
- BIO-RAD (DiaMed) ID-Card "Liss/Coombs"

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to each marked micro tube.
3. Add 25 µL of appropriate reagent to each micro tube.
4. Incubate card at +37°C for 15 minutes in an adequate incubator.
5. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
- 6 Check macroscopically for agglutination within 30 minutes.
7. Document the result.

Ortho BioVue® System (manual method)

1. Preper 3-5 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
2. Add 40 µL of appropriate reagent to the appropriate marked reaction chambers.
3. Add 10 µL of appropriate cell suspension to the reaction chambers.
4. Incubate card at +37°C for 15 minutes in an adequate incubator.
5. Centrifuge the cassette in appropriate centrifuge with the for this centrifuge unchangeable g-force.
6. The test results should be read directly after the end of centrifugation
7. Document the result.

Ortho BioVue® System (Ortho Vision® Analyzer)

The use of the mentioned above reagents is validated on the Ortho Vision® Analyzer

The Ortho Vision® Analyzer offers the following options / settings:

- 50 µL 0,8% suspension of red blood cells in PBS buffer or
- 10 µL 3-5% suspension of red blood cells in PBS buffer

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Carry out the reading and interpretation of the results according to the instruction for use of the corresponding card.

The interpretation of the results from Ortho Vision® Analyzer is performed automatically by the software.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression on human red blood cells, reactivity of the reagents mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or Variants of antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity, as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
8. Red blood cells with a positive direct Coombs-test may cause false-positive reactions.
9. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵
10. Information on limits in the instructions for use of cards used must be observed.
11. When used on the Ortho Vision® Analyzer, the test serum may remain uninterruptedly on the device (+15 to +30 °C) for a maximum of 7 days. On the other hand, if the test serum is stored overnight at +2 to +8 °C for 12 hours, the period will be extended accordingly to 14 days.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the products was conducted.

Different samples material was used and compared with a reference method.

Anti-Jk^a agreement of the results and performance data.

When testing with the Ortho monospecific IgG coombs card and the Ortho polyspecific coombs card Blood bags with the additive solutions PAGGS-M, as well as EDTA- und citrate-blood were used.

When testing with the Grifols Coombs card and the BioRad Liss/Coombs card Blood bags with the additive solutions SAG-M and PAGGS-M, as well as EDTA- und citrat blood were used.

Method	Positive	Negative	Sensitivity	Specificity
Grifols Coombs card	245 / 245	93 / 94	100%	98,94%
Bio-Rad (DiaMed) Liss/Coombs card	81 / 81	27 / 27	100%	100%
Ortho monospezifische IgG Coombs card polyspezifisch Coombs card	54 / 54	20 / 20	100%	100%

Anti-Jk^b agreement of the results and performance data:

When testing Blood bags with the additive solutions SAG-M and PAGGS-M, as well as EDTA- und citrat blood were used.

Method	Positive	Negative	Sensitivity	Specificity
Grifols Coombs card	54 / 54	20 / 20	100%	100%
Ortho monospezifische IgG Coombs card polyspezifisch Coombs card	54 / 54	20 / 20	100%	100%

LITERATUR

1. Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
2. CLSI, I/IA33-A Validatin of Automated System for Immunhematological Testing Before Implementation; Approved Guidline
CLINICAL AND LABORATORY STANDARDS INSTITUTE December 2009
3. Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. edition, Montgomery Scientific Publications 1998
4. Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3.edition, Springer-Verlag 2004.
5. Blood Transfusion Management for Patients Treated with Anti-CD38 monoclonal Antibodies.
Frontiers in Immunology November 2018 / Articiel 2616

REF Product Code	LOT Batch
 Store from - to	 Expiration Date
IVD In-Vitro Diagnostic	 EU CE symbol 0483
 Manufacturer according to 98/79/EG	 Consult instructions for use

REF

02.184-02	Anti-Jk ^a for card method Coombs-reactive, polyclonal, human	2 ml
02.184-03	Anti-Jk ^b for card method Coombs-reactive, polyclonal, human	3 ml
02.184-05	Anti-Jk ³ for card method Coombs-reactive, polyclonal, human	5 ml
02.185-02	Anti-Jk ^b for card method Coombs-reactive, polyclonal, human	2 ml
02.185-03	Anti-Jk ^a for card method Coombs-reactive, polyclonal, human	3 ml
02.185-05	Anti-Jk ³ for card method Coombs-reactive, polyclonal, human	5 ml

CE 0483

 ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Germany
730-15-5207 Version 007 / 15. January 2021