

WE LIVE TO EXCITE



ESTABLISHED
AND
EFFECTIVE

SOME THINGS SHOULD NEVER CHANGE
SEROLOGICAL HLA TYPING AND
ANTIBODY DIAGNOSTICS

HISTO TRAY
SeraScreen Abs

BAGO
DIAGNOSTICS

HISTO TRAY – THE IDEAL SOLUTION TO EXCLUDE NULL-ALLELES

Micro test trays with pre-dropped anti sera and controls for typing of HLA-ABC

- The microlymphocytotoxicity test (MLCT) is a secure and reliable confirmational test for null-alleles
- Due to ready-to-use typing trays (incl. rabbit complement) the MLCT is a simple diagnostic method
- HISTO TRAY typing trays are cost-effective for rapidly required typing results

SeraScreen Abs – EASY AND RELIABLE HLA ANTIBODY DETECTION

Micro test trays for antibody screening and detection containing frozen lymphocytes typed for HLA-ABC

- Reliable antibody diagnostic for preparation and clinical monitoring of solid organ transplantation
- Ready-to-use micro test trays with pre-dropped peripheral lymphocytes from different donors (incl. rabbit complement)

**KEEP THE
PROVEN
SeraScreen Abs
AVAILABLE
AGAIN!**

PRODUCTS	PACKAGE	REF
HISTO TRAY – Micro test trays with pre-dropped anti sera and controls for HLA typing		
HISTO TRAY AB 120 (5) (incl. 5 x 1 ml rabbit complement) for typing one specimen with two trays	5 x 2 panels	7042
HISTO TRAY ABC 72 (10) (incl. 5 x 1 ml rabbit complement) for typing one specimen per tray	10 panels	7022
HISTO TRAY ABC 120 (5) (incl. 5 x 1 ml rabbit complement) for typing one specimen with two trays	5 x 2 panels	7013
HISTO TRAY ABC 144 (5) (incl. 5 x 1 ml rabbit complement) for typing one specimen with two trays	5 x 2 panels	7035
SeraScreen Abs – Frozen cell tray HLA-class I for HLA antibody screening and identification		
SeraScreen Abs (incl. rabbit complement) for typing one specimen per tray	6 panels	7252
Reagents for HLA Serology		
HISTOPREP Lymphocyte separation medium	1 x 250 ml	70125
RPMI 1640 Cell culture medium	1 x 250 ml	70126
TERASAKI Park MEDIUM (+ 0,5% Fetal Bovine Serum) Cell culture medium	1 x 100 ml	7028
Rabbit Complement RUO (lyophilised)	10 x 1 ml	7019
Accessories		
Anti-HLA Positive Control RUO	2 ml	61100
Anti-HLA Negative Control RUO	2 ml	61101

WE LIVE TO EXCITE

EXCELLENCE, PARTNERSHIP and **PASSION** have always been the basis of our actions. Moreover, since 1947 our wish to inspire you by understanding your needs has been absolutely central – and therefore, we only develop diagnostics with added value and significance. We combine decades of experience and a deeply rooted exploratory urge for the development of our innovative products.

Thanks to our independence and global orientation, we can conduct targeted research for the development of top performance products and services in the areas of **TRANSPLANTATION, TRANSFUSION, HUMAN GENETICS** and **INFECTIOUS DISEASES.**

Our motivation is to exceed the expectations of our customers and partners. That is why we support you with exceptional service and new ways of thinking – competently, quickly and precisely.

BAG Diagnostics GmbH

Amtsgerichtsstr. 1-5
35423 Lich/Germany

T. +49 (0) 6404/925-0
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

A BAG Group company

Customer Service

T. +49 (0) 6404/925-125
F. +49 (0) 6404/925-421
M. service@bag-diagnostics.com

BAG 
DIAGNOSTICS